



# **BACTYNEA SPRAY**



## **BENEFITS**

- Can be used when infants are present
- Safe for surfaces, even the most fragile
- Bactericidal, yeasticidal et virucidal
- Foodstuff contact approved







#### **APPLICATION**

Cleaning and disinfection of non-invasive medical devices (equipment, etc.) by contact.

Researched especially for neonatal services and, more specifically, for disinfecting sensitive surfaces (sterile bubbles, insulators, incubators, etc.) when infants are present.

IIA CLASS.

#### **INSTRUCTIONS FOR USE**

Spray BACTYNEA PRET A L'EMPLOI. The foam fades quickly to facilitate wiping. Let dry. Do not rinse except in case of contact with surfaces that may come in contact with foodstuffs. Does not replace disinfection by soaking.

Complies with the decree of 19/12/13 amending the decree of 08/09/1999 concerning products for cleaning materials and articles intended to come into contact with foodstuffs.

## **CHARACTERISTICS**

Appearance: liquid
Colour: colourless
Perfume: perfume free

Density: = 1 low base Expiration date:  $\square$  2 year

Storage:



#### MICROBIOLOGICAL ACTIVITY

**Bactericidal**: EN 13727 (dirtiness conditions – 2 min – 20°C)

**Levuricidal**: EN 13624 (dirtiness conditions – 2 min – 20°C)

**Sporicidal**: EN 13704 (cleanliness conditions – 60 min – 20°C)

**Virucidal**: EN 14476 sur PRV, BVDV, H1N1, Herpès, Vaccinia virus, Rotavirus (dirtiness conditions – 2 min – 20°C), Norovirus (dirtiness conditions – 15 min – 20°C) \* Swine Pseudorabies Virus (PRV) and Bovine Viral Diarrhea Virus (BVDV) are surrogate of HBV and HCV.

## **PACKAGING**

750ml spray Réf. LB920280

## **SECURITY**

For professional use only.
For more information please refer to the MSDS.





